

KD-Validated Anti-AIFM1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1299

Product Overview

Name	KD-Validated Anti-AIFM1 Rabbit Monoclonal Antibody
Accession(Primary)	O95831
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human AIF
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 67 kDa; observed, 67 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000Flow Cytometry (FCM): 1:2,000Immunocytochemistry (IC): 1:100-1:1,000


Target information (O95831)

Gene Name	AIFM1
About Gene Symbol	AIFM1 (Apoptosis Inducing Factor Mitochondria Associated 1) is a metabolic enzyme involved in apoptosis. It is an important research target in hearing loss and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-AIFM1 antibody (Cat#ABG1299). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-AIFM1 antibody (Cat#ABG1299, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-AIFM1 antibody (Cat#ABG1299). AIFM1 expression in wild-type (WT) and AIFM1 knockdown (KD) HeLa cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-AIFM1 antibody (Cat#ABG1299, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of AIFM1 expression in HepG2 cells using anti-AIFM1 antibody (Cat#ABG1299, 1:2,000). Green, isotype control; red, AIFM1.



Immunocytochemical staining of HepG2 cells with anti-AIFM1 antibody (Cat#ABG1299, 1:1,000).

Nuclei were stained blue with DAPI; AIFM1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar, 20 μ m.