

KD-Validated Anti-Enolase 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1301

Product Overview

Name	KD-Validated Anti-Enolase 1 Rabbit Monoclonal Antibody
Accession(Primary)	P06733
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human ENO1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 47 kDa; observed, 47 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000Flow Cytometry (FCM): 1:2,000Immunocytochemistry (IC): 1:100-1:1,000


Target information (P06733)

Gene Name	ENO1
About Gene Symbol	ENO1 (Enolase 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-enolase 1 antibody (Cat#ABG1301). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-enolase 1 antibody (Cat#ABG1301, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-enolase 1 antibody (Cat#ABG1301). Enolase 1 expression in wild-type (WT) and enolase 1 (ENO1) knockdown (KD) 293T cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-enolase 1 antibody (Cat#ABG1301, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of Enolase 1 expression in HepG2 cells using anti-Enolase 1 antibody (Cat#ABG1301, 1:2,000). Green, isotype control; red, Enolase 1.



Immunocytochemical staining of HepG2 cells with anti-Enolase 1 antibody (Cat#ABG1301, 1:1,000). Nuclei were stained blue with DAPI; Enolase 1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.