

KD-Validated Anti-AGL Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1305

Product Overview

Name	KD-Validated Anti-AGL Rabbit Monoclonal Antibody
Accession(Primary)	P35573
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human AGL
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 175 kDa, observed, 150 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P35573)

Gene Name	AGL
About Gene Symbol	AGL (Amylo-Alpha-1, 6-Glucosidase, 4-Alpha-Glucanotransferase) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellula