

KD-Validated Anti-FLI1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1309

Product Overview


Name	KD-Validated Anti-FLI1 Rabbit Monoclonal Antibody
Accession(Primary)	Q01543
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human FLI1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 51 kDa; observed, 51 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000Flow Cytometry (FCM): 1:2,000Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q01543)


Gene Name	FLI1
About Gene Symbol	FLI1 (Fli-1 Proto-Oncogene, ETS Transcription Factor) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC



Immunocytochemical staining of HT-1080 cells with anti-FLI1 antibody (Cat#ABG1309, 1:1,000). Nuclei were stained blue with DAPI; FLI1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.