

KD-Validated Anti-AMBP Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1315

Product Overview

Name	KD-Validated Anti-AMBP Rabbit Monoclonal Antibody
Accession(Primary)	P02760
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human AMBP
Species Reactivity	Human
Calculated MW	Predicted, 39 kDa , observed, 39 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (P02760)

Gene Name	AMBP
About Gene Symbol	AMBP (Alpha-1-Microglobulin/Bikunin Precursor) is a formin actin nucleator involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

