

KD-Validated Anti-Cell Division Cycle 45 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1327

Product Overview

Name	KD-Validated Anti-Cell Division Cycle 45 Rabbit Monoclonal Antibody
Accession(Primary)	O75419
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CDC45
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 66 kDa; observed, 65 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000


Target information (O75419)

Gene Name	CDC45
About Gene Symbol	CDC45 (Cell Division Cycle 45) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-cell division cycle 45 antibody (Cat#ABG1327). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-cell division cycle 45 antibody (Cat#ABG1327, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-cell division cycle 45 antibody (Cat#ABG1327). Cell division cycle 45 expression in wild type (WT) and cell division cycle 45 (CDC45) knockdown (KD) HSHC cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-cell division cycle 45 antibody (Cat#ABG1327, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of Cell division cycle 45 expression in HepG2 cells using anti-Cell division cycle 45 antibody (Cat#ABG1327, 1:2,000). Green, isotype control; red, Cell division cycle 45.



Immunocytochemical staining of HepG2 cells with anti-Cell division cycle 45 antibody (Cat#ABG1327, 1:1,000). Nuclei were stained blue with DAPI; Cell division cycle 45 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.