

KD-Validated Anti-Nuclear Factor Kappa B Subunit 2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1336

Product Overview


Name	KD-Validated Anti-Nuclear Factor Kappa B Subunit 2 Rabbit Monoclonal Antibody
Accession(Primary)	Q00653
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human NFkB p100 / p52
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 97 kDa; observed, 120 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q00653)


Gene Name	NFKB2
About Gene Symbol	NFKB2 (Nuclear Factor Kappa B Subunit 2) is a transcription factor involved in NF-κB signaling and adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic ta



Copyright ©2025 Genuin Biotechnologies LLC




Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC

Nuclear factor kappa B subunit 2-Alexa Fluor® 647



Immunocytochemical staining of HepG2 cells with Nuclear factor kappa B subunit 2 antibody (Cat#ABG1336, 1:1,000). Nuclei were stained blue with DAPI; Nuclear factor kappa B subunit 2 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.