

KD-Validated Anti-DEK Rabbit Monoclonal Antibody

Rabbit monoclonal antibody





ABG1338


Product Overview

Name	KD-Validated Anti-DEK Rabbit Monoclonal Antibody
Accession(Primary)	P35659
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human DEK
Species Reactivity	Human, Rat
Calculated MW	Predicted, 43 kDa , observed, 43 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P35659)

Gene Name	DEK
About Gene Symbol	DEK (DEK Proto-Oncogene) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HeLa cells with DEK antibody (Cat#ABG1338, 1:1,000). Nuclei were stained blue with DAPI; DEK was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 μ m.