

KD-Validated Anti-Phospho-CDC37 (S13) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1339

Product Overview


Name	KD-Validated Anti-Phospho-CDC37 (S13) Rabbit Monoclonal Antibody
Accession(Primary)	Q16543
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-CDC37 (S13)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 44 kDa, observed, 44 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (Q16543)


Gene Name	CDC37
About Gene Symbol	CDC37 (Cell Division Cycle 37, HSP90 Cochaperone) is a molecular chaperone involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellu




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HT-1080 cells with Phospho-CDC37 (S13) antibody (Cat#ABG1339, 1:1,000). Nuclei were stained blue with DAPI; Phospho-CDC37 (S13) was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.