

KD-Validated Anti-Phospho-Cdk2 (Y15) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1340

Product Overview


Name	KD-Validated Anti-Phospho-Cdk2 (Y15) Rabbit Monoclonal Antibody
Accession(Primary)	P24941
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide of human Phospho-Cdk2 (Y15)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 34 kDa , observed, 34 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P24941)


Gene Name	CDK2
About Gene Symbol	CDK2 (Cyclin Dependent Kinase 2) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC