

KD-Validated Anti-Calsyntenin 1 Rabbit Monoclonal

Rabbit monoclonal antibody

ABG1350

Product Overview

Name	KD-Validated Anti-Calsyntenin 1 Rabbit Monoclonal Antibody
Accession(Primary)	O94985
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CLSTN1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 110 kDa , observed, 150 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (O94985)

Gene Name	CLSTN1
About Gene Symbol	CLSTN1 (Calsyntenin 1) is a cell junction/adhesion protein involved in cell signaling and proliferation. It is an important research target in neurodegenerative disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

