

KD-Validated Anti-Phospho-BCAR1 (Y410) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1363

Product Overview


Name	KD-Validated Anti-Phospho-BCAR1 (Y410) Rabbit Monoclonal Antibody
Accession(Primary)	P56945
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-BCAR1 (Y410)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 93 kDa , observed, 130 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P56945)


Gene Name	BCAR1
About Gene Symbol	BCAR1 (BCAR1 Scaffold Protein, Cas Family Member) is a signaling scaffold protein involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC