

KD-Validated Anti-Phospho-EGFR (Y1173) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1365

Product Overview

Name	KD-Validated Anti-Phospho-EGFR (Y1173) Rabbit Monoclonal Antibody
Accession(Primary)	P00533
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-EGFR (Y1173)
Species Reactivity	Human, Rat
Calculated MW	Predicted, 134 kDa , observed, 175 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000


Target information (P00533)

Gene Name	EGFR
About Gene Symbol	EGFR (Epidermal Growth Factor Receptor) is a receptor tyrosine kinase involved in immune signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Phospho-EGFR (Y1173) antibody (Cat#ABG1365). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Phospho-EGFR (Y1173) antibody (Cat#ABG1365, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-Phospho-EGFR (Y1173) antibody (Cat#ABG1365). Phospho-EGFR (Y1173) expression in wild type (WT) and EGFR shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-Phospho-EGFR (Y1173) antibody (Cat#ABG1365, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of Phospho-EGFR (Y1173) expression in HepG2 cells using Phospho-EGFR (Y1173) antibody (Cat#ABG1365, 1:2,000). Green, isotype control; red, Phospho-EGFR (Y1173).



Immunocytochemical staining of HepG2 cells with anti-Phospho-EGFR (Y1173) antibody (Cat #ABG1365, 1:1,000). Nuclei were stained blue with DAPI; Phospho-EGFR (Y1173) was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.