

KD-Validated Anti-Phospho-eIF4B (S406) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1366

Product Overview


Name	KD-Validated Anti-Phospho-eIF4B (S406) Rabbit Monoclonal Antibody
Accession(Primary)	P23588
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-eIF4B (S406)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 69 kDa , observed, 75 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P23588)


Gene Name	eIF4B
About Gene Symbol	eIF4B (Eukaryotic Translation Initiation Factor 4B) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of Phospho-eIF4B (S406) expression in C2C12 cells using Phospho-eIF4B (S406) antibody (ABG1366, 1:2,000). Green, isotype control; red, Phospho-eIF4B (S406).