

KD-Validated Anti-RAB5C Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1374

Product Overview

Name	KD-Validated Anti-RAB5C Rabbit Monoclonal Antibody
Accession(Primary)	P51148
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human RAB5C / RABL
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 23 kDa; observed, 23 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000Flow Cytometry (FCM): 1:2,000Immunocytochemistry (IC): 1:100-1:1,000


Target information (P51148)

Gene Name	RAB5C
About Gene Symbol	RAB5C (RAB5C, Member RAS Oncogene Family) is a small GTPase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-RAB5C antibody (Cat#ABG1374). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-RAB5C antibody (Cat#ABG1374, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-RAB5C antibody (Cat#ABG1374). RAB5C expression in wild type (WT) and RAB5C knockdown (KD) HSHC cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-RAB5C antibody (Cat#ABG1374, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of RAB5C expression in HepG2 cells using anti-RAB5C antibody (Cat#ABG1374, 1:2,000). Green, isotype control; red, RAB5C.



Immunocytochemical staining of HepG2 cells with anti-RAB5C antibody (Cat#ABG1374, 1:1,000).

Nuclei were stained blue with DAPI; RAB5C was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.