

KD-Validated Anti-Phospho-eIF4E (S209) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1376

Product Overview


Name	KD-Validated Anti-Phospho-eIF4E (S209) Rabbit Monoclonal Antibody
Accession(Primary)	P06730
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-eIF4E (S209)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 25 kDa; observed, 25 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:2,000-1:10,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P06730)


Gene Name	EIF4E
About Gene Symbol	EIF4E (Eukaryotic Translation Initiation Factor 4E) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-phospho-eIF4E (S209) antibody (Cat#ABG1376). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-phospho-eIF4E (S209) antibody (Cat#ABG1376, 1:10,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.

Western blotting analysis using anti-phospho-eIF4E (S209) antibody (Cat#ABG1376). Phospho-eIF4E (S209) expression in wild type (WT) and EIF4E knockdown (KD) HSHC cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-phospho-eIF4E (S209) antibody (Cat#ABG1376, 1:10,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.

Flow cytometric analysis of Phospho-eIF4E (S209) expression in C2C12 cells using anti-Phospho-eIF4E (S209) antibody (Cat#ABG1376, 1:2,000). Green, isotype control; red, Phospho-eIF4E (S209).



Immunocytochemical staining of C2C12 cells with Phospho-eIF4E (S209) antibody (Cat#ABG1376, 1:1,000). Nuclei were stained blue with DAPI; Phospho-eIF4E (S209) was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar, 20 μ m.