

KD-Validated Anti-Phospho-Src (Y419) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1391

Product Overview


Name	KD-Validated Anti-Phospho-Src (Y419) Rabbit Monoclonal Antibody
Accession(Primary)	P12931
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-Src (Y419)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 60 kDa , observed, 60 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (P12931)


Gene Name	SRC
About Gene Symbol	SRC (SRC Proto-Oncogene, Non-Receptor Tyrosine Kinase) is a receptor tyrosine kinase involved in receptor-mediated signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC