

KD-Validated Anti-Phospho-STAT1 (S727) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1409

Product Overview


Name	KD-Validated Anti-Phospho-STAT1 (S727) Rabbit Monoclonal Antibody
Accession(Primary)	P42224
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-STAT1 (S727)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 87 kDa , observed, 87 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P42224)


Gene Name	STAT1
About Gene Symbol	STAT1 (Signal Transducer And Activator Of Transcription 1) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cell




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with Phospho-STAT1 (S727) antibody (Cat#ABG1409, 1:1,000). Nuclei were stained blue with DAPI; Phospho-STAT1 (S727) was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 μ m.