

KD-Validated Anti-Phospho-STAT3 (Tyr705) Rabbit Monoclonal Antibody


ABG1411

Product Overview


Name	KD-Validated Anti-Phospho-STAT3 (Tyr705) Rabbit Monoclonal Antibody
Accession(Primary)	P40763
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-STAT3 (Tyr705)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 88 kDa , observed, 88 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (P40763)


Gene Name	STAT3
About Gene Symbol	STAT3 (Signal Transducer And Activator Of Transcription 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC

Phospho-STAT3 (Tyr705)-Alexa Fluor® 647