

KD-Validated Anti-Ubiquitin specific peptidase 13 Rabbit Monoclonal Antibody


ABG1417

Product Overview

Name	KD-Validated Anti-Ubiquitin specific peptidase 13 Rabbit Monoclonal Antibody
Accession(Primary)	Q92995
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human USP13
Species Reactivity	Human, Rat
Calculated MW	Predicted, 97 kDa , observed, 97 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000


Target information (Q92995)

Gene Name	USP13
About Gene Symbol	USP13 (Ubiquitin Specific Peptidase 13) is a deubiquitinase involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Ubiquitin specific peptidase 13 antibody (Cat#ABG1417). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Ubiquitin specific peptidase 13 antibody (Cat#ABG1417, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.




Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Ubiquitin specific peptidase 13 antibody (Cat#ABG1417). Ubiquitin specific peptidase 13 expression in wild type (WT) and Ubiquitin specific peptidase 13 shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-Ubiquitin specific peptidase 13 antibody (Cat#ABG1417, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of Ubiquitin specific peptidase 13 expression in HepG2 cells using Ubiquitin specific peptidase 13 antibody (Cat#ABG1417, 1:2,000). Green, isotype control; red, Ubiquitin specific peptidase 13.



Immunocytochemical staining of HepG2 cells with Ubiquitin specific peptidase 13 antibody (Cat#ABG1417, 1:1,000). Nuclei were stained blue with DAPI; Ubiquitin specific peptidase 13 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.