

KD-Validated Anti-VPS35 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1419


Product Overview

Name	KD-Validated Anti-VPS35 Rabbit Monoclonal Antibody
Accession(Primary)	Q96QK1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human VPS35
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 92 kDa, observed, 77 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (Q96QK1)

Gene Name	VPS35
About Gene Symbol	VPS35 (VPS35 Retromer Complex Component) is a protein kinase and signaling regulator involved in vesicular trafficking. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of C2C12 cells with VPS35 antibody (Cat#ABG1419, 1:1,000). Nuclei were stained blue with DAPI; VPS35 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.