

KD-Validated Anti-EWS RNA Binding Protein 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1426

Product Overview

Name	KD-Validated Anti-EWS RNA Binding Protein 1 Rabbit Monoclonal Antibody
Accession(Primary)	Q01844
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human EWSR1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 68 kDa; observed, 85 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000Flow Cytometry (FCM): 1:2,000


Target information (Q01844)

Gene Name	EWSR1
About Gene Symbol	EWSR1 (EWS RNA Binding Protein 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-EWS RNA binding protein 1 antibody (Cat#ABG1426). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-EWS RNA binding protein 1 antibody (Cat#ABG1426, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-EWS RNA binding protein 1 antibody (Cat#ABG1426). EWS RNA binding protein 1 expression in wild-type (WT) and EWS RNA binding protein 1 (EWSR1) knockdown (KD) 293T cells with 20 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-EWS RNA binding protein 1 antibody (Cat#ABG1426, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of EWS RNA binding protein 1 expression in HepG2 cells using anti-EWS RNA binding protein 1 antibody (Cat#ABG1426, 1:2,000). Green, isotype control; red, EWS RNA binding protein 1.

EWS RNA binding protein 1-Alexa Fluor® 647