

KD-Validated Anti-Histone H3 (acetyl K9) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1443

Product Overview


Name	KD-Validated Anti-Histone H3 (acetyl K9) Rabbit Monoclonal Antibody
Accession(Primary)	P68431
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Histone H3 (acetyl K9)
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 15 kDa , observed, 15 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P68431)


Gene Name	H3C1
About Gene Symbol	H3C1 (H3 Clustered Histone 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to el



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC




Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Histone H3 (acetyl K9) antibody (Cat#ABG1443). Total cell lysates (5 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Histone H3 (acetyl K9) antibody (Cat#ABG1443, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.

Western blotting analysis using anti-Histone H3 (acetyl K9) antibody (Cat#ABG1443). Histone H3 (acetyl K9) expression in wild type (WT) and Histone H3 (acetyl K9) shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-Histone H3 (acetyl K9) antibody (Cat#ABG1443, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.

Flow cytometric analysis of Histone H3 (acetyl K9) expression in HAP-1 cells using Histone H3 (acetyl K9) antibody (Cat#ABG1443, 1:2,000). Green, isotype control; red, Histone H3 (acetyl K9).



Immunocytochemical staining of HAP-1 cells with Histone H3 (acetyl K9) antibody (Cat#ABG1443, 1:1,000). Nuclei were stained blue with DAPI; Histone H3 (acetyl K9) was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 μ m.