

KD-Validated Anti-High mobility group box 2 Rabbit Monoclonal Antibody

ABG1447

Product Overview

Name	KD-Validated Anti-High mobility group box 2 Rabbit Monoclonal Antibody
Accession(Primary)	P26583
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human HMGB2
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 24 kDa , observed, 24 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (P26583)

Gene Name	HMGB2
About Gene Symbol	HMGB2 (High Mobility Group Box 2) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

