

KD-Validated Anti-QKI Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1459


Product Overview

Name	KD-Validated Anti-QKI Rabbit Monoclonal Antibody
Accession(Primary)	Q96PU8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human QK1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 38 kDa , observed, 38 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q96PU8)

Gene Name	QKI
About Gene Symbol	QKI (QKI, KH Domain Containing RNA Binding) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of C2C12 cells with QKI antibody (Cat#ABG1459, 1:1,000). Nuclei were stained blue with DAPI; QKI was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.