

KD-Validated Anti-Nesprin3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1465

Product Overview

Name	KD-Validated Anti-Nesprin3 Rabbit Monoclonal Antibody
Accession(Primary)	Q6ZMZ3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Nesprin3
Species Reactivity	Human
Calculated MW	Predicted, 112 kDa, observed, 100 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (Q6ZMZ3)

Gene Name	SYNE3
About Gene Symbol	SYNE3 (Spectrin Repeat Containing Nuclear Envelope Family Member 3) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

