

KD-Validated Anti-Synaptophysin like 1 Rabbit Monoclonal

Rabbit monoclonal antibody




ABG1466


Product Overview

Name	KD-Validated Anti-Synaptophysin like 1 Rabbit Monoclonal Antibody
Accession(Primary)	Q16563
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human SYPL1
Species Reactivity	Human
Calculated MW	Predicted, 29 kDa , observed, 29-37 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (Q16563)

Gene Name	SYPL1
About Gene Symbol	SYPL1 (Synaptophysin Like 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HeLa cells with Synaptophysin like 1 antibody (Cat#ABG1466, 1:1,000). Nuclei were stained blue with DAPI; Synaptophysin like 1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.