

KD-Validated Anti-Tribbles pseudokinase 3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1467

Product Overview

Name	KD-Validated Anti-Tribbles pseudokinase 3 Rabbit Monoclonal Antibody
Accession(Primary)	Q96RU7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TRIB3
Species Reactivity	Human
Calculated MW	Predicted, 40 kDa , observed, 40 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (Q96RU7)

Gene Name	TRIB3
About Gene Symbol	TRIB3 (Tribbles Pseudokinase 3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

