

KD-Validated Anti-CRKL Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1473

Product Overview


Name	KD-Validated Anti-CRKL Rabbit Monoclonal Antibody
Accession(Primary)	P46109
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CrkL
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 34 kDa , observed, 37 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P46109)

Gene Name	CRKL
About Gene Symbol	CRKL (CRK Like Proto-Oncogene, Adaptor Protein) is a signaling adaptor protein involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blotting analysis using anti-CRKL antibody (Cat#ABG1473). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-CRKL antibody (Cat#ABG1473, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-CRKL antibody (Cat#ABG1473). CRKL expression in wild type (WT) and CRKL shRNA knockdown (KD) HT-1080 cells with 30 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-CRKL antibody (Cat#ABG1473, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.