

KD-Validated Anti-Phospho-PKC alpha (T497) Rabbit Monoclonal Antibody


ABG1476

Product Overview


Name	KD-Validated Anti-Phospho-PKC alpha (T497) Rabbit Monoclonal Antibody
Accession(Primary)	P17252
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-PKC alpha (T497)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 77 kDa, observed, 73 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (P17252)


Gene Name	PRKCA
About Gene Symbol	PRKCA (Protein Kinase C Alpha) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in aging research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC