

KD-Validated Anti-RUNX1/RUNX2/RUNX3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1477

Product Overview


Name	KD-Validated Anti-RUNX1/RUNX2/RUNX3 Rabbit Monoclonal Antibody
Accession(Primary)	Q01196/Q13761/Q13950
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human RUNX1/RUNX2/RUNX3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 49 kDa ; Observed, 51 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Flow cytometry (FCM): 1:2,000

Target information (Q01196/Q13761/Q13950)


Gene Name	RUNX1/RUNX2/RUNX3
About Gene Symbol	RUNX1/RUNX2/RUNX3 (RUNX Family Transcription Factor 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanism.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC