

KD-Validated Anti-Histone deacetylase 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1478

Product Overview


Name	KD-Validated Anti-Histone deacetylase 1 Rabbit Monoclonal Antibody
Accession(Primary)	Q13547
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human HDAC1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 55 kD a , observed, 55 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q13547)


Gene Name	HDAC1
About Gene Symbol	HDAC1 (Histone Deacetylase) is a histone deacetylase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024, Genuin Biotechnologies LLC



Copyright ©2024, Genuin Biotechnologies LLC



Copyright ©2024, Genuin Biotechnologies LLC