

KD-Validated Anti-HNRNPK Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1479

Product Overview


Name	KD-Validated Anti-HNRNPK Rabbit Monoclonal Antibody
Accession(Primary)	P61978
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human hnRNP K
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 51 kDa, observed, 60 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:10,000 Flow Cytometry (FCM): 1:2,000

Target information (P61978)


Gene Name	HNRNPK
About Gene Symbol	HNRNPK (Heterogeneous Nuclear Ribonucleoprotein K) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC