

KD-Validated Anti-Katanin catalytic subunit A1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1481

Product Overview


Name	KD-Validated Anti-Katanin catalytic subunit A1 Rabbit Monoclonal Antibody
Accession(Primary)	O75449
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human p60 katanin
Species Reactivity	Human
Calculated MW	Predicted, 56 kDa, observed, 49 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (O75449)


Gene Name	KATNA1
About Gene Symbol	KATNA1 (Katanin Catalytic Subunit A1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC