

KD-Validated Anti-Galectin 3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1482


Product Overview

Name	KD-Validated Anti-Galectin 3 Rabbit Monoclonal Antibody
Accession(Primary)	P17931
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Galectin 3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 26 kDa, observed, 26 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P17931)

Gene Name	LGALS3
About Gene Symbol	LGALS3 (Galectin 3) is an extracellular matrix laminin protein involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with anti-Galectin 3 antibody (Cat#ABG1482, 1:1,000). Nuclei were stained blue with anti-DAPI; Galectin 3 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.