

KD-Validated Anti-Phospho-p38 (T180) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1484

Product Overview


Name	KD-Validated Anti-Phospho-p38 (T180) Rabbit Monoclonal Antibody
Accession(Primary)	Q16539
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-p38 (T180)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 41 kD a , observed, 38 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q16539)


Gene Name	MAPK14
About Gene Symbol	MAPK14 (Mitogen-Activated Protein Kinase 14) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and JAK-STAT/cytokine signaling. It is an important research target in inflammatory disease and is widely studied to elucidate cellular mecha



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC