

KD-Validated Anti-gamma Tubulin Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1490

Product Overview

Name	KD-Validated Anti-gamma Tubulin Rabbit Monoclonal Antibody
Accession(Primary)	P23258
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Tubulin gamma
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 51 kDa, observed, 48 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P23258)

Gene Name	TUBG1
About Gene Symbol	TUBG1 (Tubulin Gamma 1) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

