

KD-Validated Anti-Thioredoxin reductase 2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1491

Product Overview


Name	KD-Validated Anti-Thioredoxin reductase 2 Rabbit Monoclonal Antibody
Accession(Primary)	Q9NNW7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TXNRD2
Species Reactivity	Human
Calculated MW	Predicted, 57 kDa, observed, 57 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry(FCM):1:2,000

Target information (Q9NNW7)


Gene Name	TXNRD2
About Gene Symbol	TXNRD2 (Thioredoxin Reductase 2) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC