

KD-Validated Anti-CCDC98 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1495

Product Overview


Name	KD-Validated Anti-CCDC98 Rabbit Monoclonal Antibody
Accession(Primary)	Q6UWZ7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CCDC98
Species Reactivity	Human
Calculated MW	Predicted, 47 kDa , observed, 50 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q6UWZ7)


Gene Name	ABRAXAS1
About Gene Symbol	ABRAXAS1 (Antibodyraxas 1, BRCA1 A Complex Subunit 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blotting analysis using anti-CCDC98 antibody (Cat#ABG1495). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-CCDC98 antibody (Cat#ABG1495, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-CCDC98 antibody (Cat#ABG1495). CCDC98 expression in wild type (WT) and CCDC98 shRNA knockdown (KD) HT-1080 cells with 30 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-CCDC98 antibody (Cat#ABG1495, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Flow cytometric analysis of CCDC98 expression in HT-1080 cells using CCDC98 antibody (Cat#ABG1495, 1:2,000). Green, isotype control; red, CCDC98.