

KD-Validated Anti-Beta-secretase 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1497

Product Overview


Name	KD-Validated Anti-Beta-secretase 1 Rabbit Monoclonal Antibody
Accession(Primary)	P56817
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human BACE1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 56 kD a , observed, 70 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (P56817)


Gene Name	BACE1
About Gene Symbol	BACE1 (Beta-Secretase 1) is a protease enzyme involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC