

KD-Validated Anti-C-terminal binding protein 2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1500

Product Overview


Name	KD-Validated Anti-C-terminal binding protein 2 Rabbit Monoclonal Antibody
Accession(Primary)	P56545
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CTBP2
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 49 kDa , observed, 47 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P56545)


Gene Name	CTBP2
About Gene Symbol	CTBP2 (C-Terminal Binding Protein 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC