

KD-Validated Anti-Diacylglycerol O-Acyltransferase 1 Rabbit Monoclonal Antibody


ABG1503

Product Overview

Name	KD-Validated Anti-Diacylglycerol O-acyltransferase 1 Rabbit Monoclonal Antibody
Accession(Primary)	O75907
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human DGAT1
Species Reactivity	Human
Calculated MW	Predicted, 55 kDa , observed, 48 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000


Target information (O75907)

Gene Name	DGAT1
About Gene Symbol	DGAT1 (Diacylglycerol O-Acyltransferase 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Diacylglycerol O-acyltransferase 1 antibody (Cat#ABG1503). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Diacylglycerol O-acyltransferase 1 antibody (Cat#ABG1503, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Diacylglycerol O-acyltransferase 1 antibody (Cat#ABG1503). Diacylglycerol O-acyltransferase 1 expression in wild type (WT) and Diacylglycerol O-acyltransferase 1 shRNA knockdown (KD) HT-1080 cells with 30 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-Diacylglycerol O-acyltransferase 1 antibody (Cat#ABG1503, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of Diacylglycerol O-acyltransferase 1 expression in HeLa cells using Diacylglycerol O-acyltransferase 1 antibody (Cat#ABG1503, 1:2,000). Green, isotype control; red, Diacylglycerol O-acyltransferase 1.