

KD-Validated Anti-Desmoglein 2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1508

Product Overview


Name	KD-Validated Anti-Desmoglein 2 Rabbit Monoclonal Antibody
Accession(Primary)	Q14126
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Desmoglein 2
Species Reactivity	Human
Calculated MW	Predicted, 122 kDa , observed, 160 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q14126)


Gene Name	DSG2
About Gene Symbol	DSG2 (Desmoglein 2) is a cell junction/adhesion protein involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC