

KD-Validated Anti-GCLM Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1522

Product Overview

Name	KD-Validated Anti-GCLM Rabbit Monoclonal Antibody
Accession(Primary)	P48507
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GCLM
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 31 kDa , observed, 28 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (P48507)

Gene Name	GCLM
About Gene Symbol	GCLM (Glutamate-Cysteine Ligase Modifier Subunit) is a biosynthetic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms

