

KD-Validated Anti-Golgi integral membrane protein 4 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1526

Product Overview

Name	KD-Validated Anti-Golgi integral membrane protein 4 Rabbit Monoclonal Antibody
Accession(Primary)	O00461
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GOLPH4
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 82 kDa , observed, 130 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (O00461)

Gene Name	GOLIM4
About Gene Symbol	GOLIM4 (Golgi Integral Membrane Protein 4) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

