

KD-Validated Anti-Golgi reassembly stacking protein 1

Rabbit monoclonal antibody


ABG1527

Product Overview

Name	KD-Validated Anti-Golgi reassembly stacking protein 1 Rabbit Monoclonal Antibody
Accession(Primary)	Q9BQQ3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GRASP65
Species Reactivity	Human
Calculated MW	Predicted, 46 kDa , observed, 65 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000


Target information (Q9BQQ3)

Gene Name	GORASP1
About Gene Symbol	GORASP1 (Golgi Reassembly Stacking Protein 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Golgi reassembly stacking protein 1 antibody (Cat#ABG1527). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Golgi reassembly stacking protein 1 antibody (Cat#ABG1527, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Golgi reassembly stacking protein 1 antibody (Cat#ABG1527). Golgi reassembly stacking protein 1 expression in wild type (WT) and Golgi reassembly stacking protein 1 shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-Golgi reassembly stacking protein 1 antibody (Cat#ABG1527, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of Golgi reassembly stacking protein 1 expression in HeLa cells using Golgi reassembly stacking protein 1 antibody (Cat#ABG1527, 1:2,000). Green, isotype control; red, Golgi reassembly stacking protein 1.