

KD-Validated Anti-TCEB2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1532

Product Overview

Name	KD-Validated Anti-TCEB2 Rabbit Monoclonal Antibody
Accession(Primary)	Q15370
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TCEB2
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 13 kDa , observed, 16 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q15370)

Gene Name	ELOB
About Gene Symbol	ELOB (Elongin B) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

