

KD-Validated Anti-GATA2/GATA3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1538

Product Overview


Name	KD-Validated Anti-GATA2/GATA3 Rabbit Monoclonal Antibody
Accession(Primary)	P23769/P23771
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GATA2/3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 48,51 kDa ; Observed , 51 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P23769/P23771)


Gene Name	GATA2/3
About Gene Symbol	GATA2/3 (GATA Binding Protein 2) is a transcription factor involved in angiogenesis. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HAP1 cells with anti-GATA2/3 antibody (Cat #ABG1538, 1:1,000). Nuclei were stained blue with DAPI; GATA2/3 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.