

KD-Validated Anti-HuR / ELAVL1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1542

Product Overview

Name	KD-Validated Anti-HuR / ELAVL1 Rabbit Monoclonal Antibody
Accession(Primary)	Q15717
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human HuR / ELAVL1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 36 kDa, observed, 30 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q15717)

Gene Name	ELAVL1
About Gene Symbol	ELAVL1 (ELAV Like RNA Binding Protein 1) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

