

KD-Validated Anti-GTF2I Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1545

Product Overview


Name	KD-Validated Anti-GTF2I Rabbit Monoclonal Antibody
Accession(Primary)	P78347
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GTF2I
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 113 kDa , observed , 135 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000

Target information (P78347)

Gene Name	GTF2I
About Gene Symbol	GTF2I (General Transcription Factor Ili) is a receptor tyrosine kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blotting analysis using anti-GTF2I antibody (Cat#ABG1545). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-GTF2I antibody (Cat#ABG1545, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-GTF2I antibody (Cat#ABG1545). GTF2I expression in wild type (WT) and GTF2I shRNA knockdown (KD) HeLa cells with 30 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-GTF2I antibody (Cat#ABG1545, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.